












Symbols Glossary

51034 ver 02

| Symbol | Title | Reference number |
|--|---|------------------|
|  | Caution | 5.4.4 |
|  | Consult instructions for use | 5.4.3 |
|  | Do not use if package is damaged | 5.2.8 |
|  | Do not re-use | 5.4.2 |
|  | Use-by date | 5.1.4 |
|  | Batch code | 5.1.5 |
|  | Sterilized using ethylene oxide | 5.2.3 |
|  | Temperature limit | 5.3.7 |
|  | Authorized Representative in the European Community | 5.1.2 |
|  | Manufacturer | 5.1.1 |
|  | Catalog number | 5.1.6 |
| <p>* According to: ISO 15223-1:2016 (Corrected version 2017-03) Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied -- Part 1: General requirements</p> | | |