












OSSIX® Regenerative Solutions Matrix

Your Talent. Our Expertise.

The OSSIX® product family provides a palette of solutions for a wide range of clinical procedures

Procedures	OSSIX® Plus Ossifying, Cross-Linked Barrier Membrane	OSSIX® Volumax Volumizing, Ossifying Collagen Scaffold	OSSIX™ Bone Ossifying Collagen Sponge	images
Lateral Bone Augmentation	✓✓ Possible with all bone grafts	✓✓ (with/without bone graft) Subject to supporting/ non- supporting bone	✓✓✓ Subject to supporting/ non-supporting bone with OSSIX Plus	
Socket Preservation	✓✓ with/without bone graft (contained defect)	✓✓ with/without bone graft	✓✓✓ no need for membrane	
Crestal Sinus Elevation	n/a	n/a	✓✓✓ no need for membrane	
Lateral Window Sinus Elevation	✓✓✓ over lateral window	✓✓✓ over lateral window	✓	
Vertical Augmentation	✓✓ with bone graft	✓✓ with bone graft	✓✓ with membrane	
Preventing / Remediating Buccal Bone Loss	✓✓ with bone graft	✓✓✓ without bone graft	✓✓ with membrane	
2 nd Stage Implant Surgery Dehiscence	✓✓ with bone graft	✓✓✓ without bone graft	✓✓ with membrane	
2 nd Stage Implant Surgery Vertical	✓✓ with bone graft	✓✓ without bone graft (Poncho)	✓✓ with membrane	
Gingival Recession	n/a	✓	n/a	
Increase Soft Tissue Volume*	n/a	✓✓	n/a	
Mandibular Nerve Repositioning*	n/a	✓✓ with/without bone graft	n/a	

✓ Possible ✓✓ Recommended ✓✓✓ Highly recommended

* Limited clinical experience with this procedure.

Please refer to IFUs for the full information on indications and contraindications.








OSSIX[®] Regenerative Solutions Matrix

Your Talent. Our Expertise.

The OSSIX[®] product family provides a palette of solutions for a wide range of clinical procedures

datumdental

For more information on OSSIX regenerative products and activities in your region:

www.ossixdental.com     

© 2018 Datum Dental, Ltd. All rights reserved. The Datum Dental logo, GLYMATRIX and OSSIX are trademarks and/or registered trademarks of Datum Dental Ltd, Datum Biotech Ltd, their subsidiaries, and/or affiliates in the United States and/or other countries. All other company or product names are the trademarks or registered trademarks of their respective holders.

Please refer to IFUs for the full information on indications and contraindications.